


## Real-Time Cell Health Monitoring, Cell-by-Cell, in 2D and 3D in the Far-Red.


In **flow cytometry** and **fluorescence microscopy** (including **high content screening**), a reliable estimation of cell viability is important because it is central to assays for **apoptosis** and ***in vitro* toxicology**. Likewise, it is often a useful measure for sample quality and for robust phenotypic analysis of clinical samples. To better report such cell viability DRAQ7™ (Edward, 2011) was developed, based on DRAQ5™. Essentially, this water-soluble probe has identical spectral properties to DRAQ5™ and as such does not overlap with the majority of visible range fluors. Likewise, RNA binding is very weak (undetectable by flow cytometry). However, the chemical modification has rendered DRAQ7™ membrane impermeant, thus it does not cross the membrane of viable cells but rapidly enters “leaky” cells and labels nuclear DNA. Therefore, DRAQ7™ can be used as a new far-red reporter of cell viability and conversely cell membrane-permeabilization resulting from damage, apoptosis and necrosis. As such, DRAQ7™ is an ideal spectrally-shifted replacement for agents such as DAPI, propidium iodide and TOTO-3 and may offer some new information and a wider assay window. For HCS studies in drug discovery and *in vitro* toxicology DRAQ7™ can be applied as a reporter of cell-membrane permeabilization, combined with live cell-permeant DNA dyes (e.g. CyTRAK Orange™ and Hoechst dyes). DRAQ7™ can be applied as a viability reporter in mitotic index assays and in studies into cell health in response to insults. As an example, for the typical HCS cell health assay, which combines a dye for “all events” (Hoechst 33342), a mitochondrial membrane potential reporter (e.g. TMRM) and a cell viability dye (TOTO-3), the latter component can be replaced with DRAQ7™ for less spectral overlap and thus wider detection windows for the TMRM and cell viability components whilst significantly reducing total assay reagent costs.



As shown left, SU-DHL-4 cells cultured in the presence of DRAQ7™ at both the standard concentration or 3.3X excess show no significant impact on growth curves compared to untreated controls making DRAQ7™ an ideal candidate as a reporter of cell death in real-time, long-term viability and toxicity assays. DRAQ7™ has been robustly tested and exemplified in key cytometry publications (Akagi et al, 2013a; Akagi et al, 2013b; Wlodkowic et al, 2013; Smith et al, 2013) for compatibility with multi-colour experiments and uniquely for long-term, real-time analysis. It has recently been utilised in

imaging procedures to monitor cell viability in 2-D and 3-D spheroid/micro-tissue assays including a study on glioblastoma-derived stem cell lines in response to a library of chemotherapeutic agents (shown below and detailed elsewhere in a separate joint BioStatus / Imagen-Biotech white paper).

Additionally, DRAQ7™ has been applied to a variety of **long-term and real-time cell health assays** that benefit from its fundamental features. For example, DRAQ7™ has been applied to monitoring cell health in 7-day culture of pancreas tissue (Marciniak et al, 2013), for real-time nanoparticle toxicity monitoring (Ware et al, 2014), for real-time study of mitochondrially-regulated apoptosis (Liang et al, 2015) and detecting subcutaneous extracellular dsDNA in better understanding its adjuvant role in vaccination (Wang et al, 2015).



Previous methodologies have attempted to positively mark cells based on metabolic competence. DRAQ7™ labels only (damaged, dying and dead) membrane-compromised cells and, importantly, red excitation minimises risk of DNA damage when capturing time-lapse images whilst reliably monitoring viability in real-time, cell-by-cell.

DRAQ7™ offers new dimensions and opportunities for performance of high value phenotypic and *in vitro* toxicity cell-based assays in drug discovery and development, that can be applied across different platforms including flow cytometry, fluorescence microscopy and high content imaging platforms.

## Reference List

Akagi, J., Kordon, M., Zhao, H., et al. (2013a). Real-time cell viability assays using a new anthracycline derivative DRAQ7. *Cytometry Part A* **83**, 227-234.

Akagi, J., Skommer, J., Matuszek, et al. (2013b). Multivariate analysis of apoptotic markers versus cell cycle phase in living human cancer cells by microfluidic cytometry. *Proc. SPIE* **8615**.

Edward, R. (2011). DRAQ7™ - A unique far-red viability dye provides new combinations for cell health & *in vitro* tox assays. Poster presentation: Society for Biomolecular Sciences Annual Conference, Orlando, USA.

Liang, J. R., Martinez, A., Lane, J. D., et al. (2015). USP30 deubiquitylates mitochondrial Parkin substrates and restricts apoptotic cell death. *EMBO reports*, e201439820.

Marciniak, A., Selck, C., Friedrich, B., & Speier, S. (2013). Mouse pancreas tissue slice culture facilitates long-term studies of exocrine and endocrine cell physiology in situ. *PLoS One* 8(11) e78706.

Smith, P. J., Wiltshire, M., Chappell, et al. Kinetic analysis of intracellular Hoechst 33342—DNA interactions by flow cytometry: Misinterpretation of side population status? *Cytometry Part A* 83.1 (2013): 161-169.

Wang, J., Li, B., & Wu, M. X. (2015). Effective and lesion-free cutaneous influenza vaccination. *Proceedings of the National Academy of Sciences*, 112(16), 5005-5010.

Ware, M. J., Godin, B., Singh, N., et al. (2014). Analysis of the influence of cell heterogeneity on nanoparticle dose response. *ACS nano*, 8(7), 6693-6700.

Wlodkowic, D., Akagi, J., Dobrucki, J., et al. (2013). Kinetic viability assays using DRAQ7 probe. *Current Protocols in Cytometry*, 9-41.

---

If you would like to know more about DRAQ7™ or any other BioStatus product get in touch ..

email: [enquiry@biostatus.com](mailto:enquiry@biostatus.com)  
phone: +44 (0)1509 558163  
website: [www.biostatus.com](http://www.biostatus.com)

BioStatus products are for research use only.  
BioStatus products are the subject of several international patents. DRAQ5™, CyTRAK Orange™ & DRAQ7™ are trademarks of BioStatus Limited.

BioStatus Limited The BioStatus logo, consisting of three orange circles of varying sizes above the word "biostatus" in green lowercase letters.

56A Charnwood Rd, Shepshed,  
Leicestershire LE12 9NP UK

Company registered in England and Wales. No. 3079239

An ISO9001:2008 registered company



© Copyright BioStatus Ltd. September 2015