www.monobind.com

100 North Pointe Drive Lake Forest, CA 92630 TEL 949.951.2665 FAX 949.951.3539

Advisory Notice

he World Resource for Diagnostic Products

Monobind Inc.

Date Released:2016-May-6Posted By:Quality ControlEffective Date:2016-Febuary-15Product Affected:T4 Enzyme (ELISA) & Tracer (CLIA) in Total T4, Free T4, NT4, & Thyroid VASTBrief Description:New Amber Vial

To whom it may concern,

Recent product improvements to Monobind T4 Assays have shown the use of Amber Vials as a better safeguard for "<u>T4 Enzyme and Tracer</u>" and "<u>Free T4 Enzyme and Tracer</u>" against light sensitivity. Product updates are summarized below:

Reagent	Previous Packaging	New Packaging	Lots Affected
T4 Enzyme	Clear Vial	Amber Vial	2Q1B6 & Newer
T4 Tracer	Clear Vial	Amber Vial	2Q1E6 & Newer
Free T4 Enzyme	Clear Vial	Amber Vial	12Q1C6 & Newer
Free T4 Tracer	Clear Vial	Amber Vial	12Q1C6 & Newer
Strep T4 Enzyme	Clear Vial	Amber Vial	82Q1A6 & Newer
Strep T4 Tracer	Clear Vial	Amber Vial	82Q1A6 & Newer
Strep Free T4 Enzyme	Clear Vial	Amber Vial	72Q1A6 & Newer
Strep Free T4 Tracer	Clear Vial	Amber Vial	72Q1A6 & Newer
Neonatal T4 Enzyme	Clear Vial	Amber Vial	26Q1A6 & Newer

These reagents are present as kit components in the following products:

 Total T4 AccuBind ELISA 	Item #: 225-300
 Total T4 AccuLite CLIA 	Item #: 275-300
Total T4 SBS AccuBind ELISA	Item #: 8225-300
 Total T4 SBS AccuLite CLIA 	Item #: 8275-300
 Free T4 AccuBind ELISA 	Item #: 1225-300
 Free T4 AccuLite CLIA 	Item #: 1275-300
Free Thyroid Panel AccuBind V	AST ELISA Item #: 7025-300
Free Thyroid Panel AccuLite VA	ST CLIA Item #: 7075-300
Thyroid Panel AccuBind VAST I	ELISA Item #: 8025-300
Thyroid Panel AccuLite VAST C	LIA Item #: 8075-300
 Neonatal T4 AccuBind ELISA 	Item #: 2425-300
other concets of the product remain	in unchanged including the IFU. T

All other aspects of the product remain unchanged, including the IFU. The Design Output has been updated to reflect the new packaging.

As an official of the company, I verify this information to be true and correct.

Sincerely,

, Ashatola

Anthony Shatola, Quality Control Director